Product

ActiMaris® Wond gel 20 gr

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General information – ActiMaris® wound gel (20g)

1. Product Profile

ActiMaris® wound gel is made of alkaline ionized sea water with sea salt (3%), NaOCl (0.2%) and active oxygen (O_2) in the form of HOCl (0,0044%) as lithium magnesium sodium silicate.

ActiMaris® wound gel in interaction with the mucous membrane, skin and wound is nascent.

Active oxygen, sodium hypochlorite and sea salt have an eliminating effect on typical sore pathogens without damaging the tissue. The decontaminating effect against MRSA / VRE tested positive. The oxidative effect of a singlet oxygen (O_2) and hypertonic saline sea are also responsible for the neutralization of wound odors.

The basic pH of ActiMaris® wound gel supported by the provision of OH⁻ ions (as part of a naturally balanced redox system), additionally controlled the release of oxygen.

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The increased salinity of ActiMaris® wound gel in turn ensures an optimum hyperosmotic effect so that peri-wound swelling subsides quickly.

2. Fields of application

ActiMaris® wound gel is suitable:

For cleaning, humidification, decontamination, biofilm resolution, swelling and physiological debridement, to create an optimal «micro-environment» for:

- acute mechanical wounds (skin lacerations, bites, cuts, abrasions, tears and bruises) and postoperative wounds.
- chronic wounds (decubitus ulcers, arterial-venous ulcers, diabetic ulcers)
- Necrotic, malodorous wounds and ulcerating tumors, even with cavities (Hollow spaces)
- Thermal and chemical wounds (burns 1st to 3rd degree)
- entry portals of urological catheters and PEG tubes and drains
- with venous ulcers

For the local treatment of skin and mucosa in inflammatory and infectious processes. For the prevention of infection at the mucosa and skin, burns and other wound types. For this, ActiMaris® wound gel is well suited as a wound filler.

3. Instructions for use

ActiMaris® wound gel should be applied at room temperature. Before this, the wound gel can be warmed to body temperature before use. Before application of the wound gel, the wound should be cleaned with ActiMaris® wound irrigation solution (Forte or Sensitiv). The wound gel is then liberally applied on or in the wound and covered with an inert compress, gauze etc. The duration and frequency of application should be adapted to the individual wound conditions and should be carried out until all deposits and necrosis are easily removed and the wound is visually clean.

Depending on the wound condition, a bandage change can be done according to indications and phases, every 2 to 3 days. In case of non-infected wounds, necrosis and coatings can be mechanically removed depending on the condition of the wound. In cases of difficult to remove bandages, it is recommended to hydrate the wound dressings with wound irrigation solution ActiMaris® (Forte or sensitivity) for at least 5 minutes until a smooth removal of the bandages is possible without trauma to the wound surface. A primary surgical therapy is recommended with thick, firmly adhering fibrinous or necrotic wound deposits. For cleaning wounds, ActiMaris® wound gel (See also Section 8: interactions) can be combined with inert dressings or other wound dressings.

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ActiMaris® wound gel is suitable for long-term (> 24h) as well as repeated use for during dressing changes, even with irritated skin and mucous membranes. The application is painless and is generally well tolerated by allergy sufferers. Inhibition of granulation and epithelialization is not to be feared. The wound gel is dermatologically assessed as safe.

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- Slight temporary burning sensation
- If applied to nose area, slight burning and sneezing

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Infants and newborns:

Due to insufficient experience on long-term exposure, ActiMaris® wound gel should only be applied after careful medical benefit-risk assessment.

9. Interactions

ActiMaris® wound gel acts as a bio-physical compound due to the content of ROS substances ($O_2 + OCl^-$) energy and oxidizing, which is why a simultaneous combination with other local, active substances and / or interactive wound dressings is in principle unnecessary. Contact with non-precious metals should be avoided. The wound can be covered with modern wound dressings after applying ActiMaris® wound gel.

On the other hand, a good compatibility with modern wound dressings such as silver dressings has been shown in clinical trials. Likewise, there is ample clinical evidence of good tolerability when combined with inert, non-active bandages and wound dressings such as compresses, gauze, pads, sponges, gels, hydro fibers, alginates, hydrocolloids etc. Concomitant use of ActiMaris® wound irrigation solution with the ActiMaris® gel is possible and useful.

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Sea salt, water, lithium magnesium sodium silicate

Clear, colorless, basic and hypertonic gel with characteristic smell; contains 3% sea salt (Sal Maris) and 0,2% NaOCl.

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simple medical solutions

-  Lot number
-  Order number
-  Follow instructions for use
-  Use by: year and month
-  Storage (temperature)
-  Protect from direct sunlight
-  Manufacturer

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General information – ActiMaris® Forte wound rinsing solution (300ml)

1. Product profile

ActiMaris® Forte wound irrigation solution is made of sea water with sea salt (3%) and NaOCl (0.2%). ActiMaris® Forte wound irrigation solution in interaction with the mucous membrane, skin and wound is nascent. Active oxygen, sodium hypochlorite and sea salt have an eliminating effect on typical sore pathogens without damaging the tissue. The decontaminating effect against MRSA / VRE tested positive. The oxidative effect of a singlet oxygen (1O_2) and hypertonic saline are also responsible for the neutralization of wound odors. The basic pH of ActiMaris® Forte wound irrigation solution supported the provision of OH⁻ ions (as part of a naturally balanced redox system) and in addition, controlled the release of oxygen. Sea salt contains sodium chloride as its main ingredient (NaCl = 98%) as well as 2% other sea salts. The increased salinity of ActiMaris® Forte wound irrigation solution in turn ensures optimum hyperosmotic effect so that peri-wound swelling subsides quickly.

2. Fields of application

ActiMaris® Forte wound rinsing solution is suitable: For cleaning, humidification, decontamination, biofilm resolution, swelling and physiological debridement, to create an optimal «micro-environment» for:

- Acute mechanical wounds (skin lacerations, bites, cuts and abrasions, crack and contused) and postoperative wounds.
- Chronic wounds (pressure ulcers, arteriovenous ulcers, diabetic ulcers)
- Necrotic, malodorous wounds and ulcerating tumors, even with cavities (Hollow spaces)
- Thermal and chemical wounds (burns 1st to 3rd degree)
- entry portals of urological catheters and PEG tubes and drains
- with venous ulcers

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During dressing changes, loosen encrusted bandages or other encrusted wound dressings.

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For the prevention of infection of the mucosa and skin, burns and other wound types.

3. Instructions for use

ActiMaris® Forte wound rinsing solution can be applied at room temperature. In the event of strong pain, it is suggested that the necessary amount be heated to body temperature in a separate bowl before use (rinsing solution is to be heated to approximately 37°C). Before using ActiMaris® Forte wound rinsing solution, the wounds should first be rinsed and cleaned with ActiMaris® Sensitive wound rinsing solution, in order to remove dirt as well as easily removable coatings. The subsequent deep cleansing and odor reduction will be implemented with ActiMaris® Forte wound rinsing solution; the soaked sterile bandages or wound dressings are applied for this purpose (see also section 8: Interactions).

The application should preferably be in the form of the classic wet / dry phase («moist balance»).

The recommended application time of the wet phase in case of severe fibrinous or purulent, infected or strong foul-smelling wounds is → 5-10 minutes.

The duration and frequency of application should be adapted to the individual wound conditions. The use of the wound rinsing solution is not limited to time and frequency. Necrosis and coatings can be removed mechanically, depending on the wound condition. Primary surgical sanitation is advised with thick, adherent, fibrous or necrotic wound coatings.

After the improvement of the wound condition, it is advised to continue further wound cleansing with ActiMaris® Sensitive wound rinsing solution until the wound is visually clean.

ActiMaris Sensitive wound irrigation solution is recommended for patients with sensitive skin.

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ActiMaris® Forte wound irrigation solution is suitable for long-term (> 24h) as well as repeated use during dressing changes even with irritated skin and mucous. The application is painless and is generally well tolerated by allergy sufferers. Inhibition of granulation and epithelialization is not to be feared. The wound gel is dermatologically assessed as safe.

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- Slight temporary burning sensation
- Slight burning and sneezing if applied to nose area,

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9. Interactions:

ActiMaris® wound gel acts as a bio-physical compound due to the content of ROS substances (1O_2 + OCl⁻) energy and oxidizing, which is why a simultaneous combination with other local, active substances and / or interactive wound dressings is in principle unnecessary. Contact with non-precious metals should be avoided. The wound can be covered with modern wound dressings after applying ActiMaris® wound gel.

On the other hand, clinical trials have shown a good compatibility with modern wound dressings such as silver dressings. Likewise, there is ample clinical evidence of good tolerability when combined with inert, non-active bandages and wound dressings such as compresses, gauze, pads, sponges, gels, hydro fibers, alginates, hydrocolloids etc. The concomitant use of ActiMaris® wound irrigation solution with the ActiMaris® gel is possible and useful.

10. General safety instructions

ActiMaris® wound gel is suitable only for external wounds on skin or mucosa. Do not use for infusion or injection! Do not swallow!

Use only closed and undamaged containers. Keep out of reach of children!

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Sea salt, water, lithium magnesium sodium silicate

Clear, colorless, basic and hypertonic gel with characteristic smell; contains 3% sea salt (Sal Maris) and 0.2% NaOCl.

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ActiMaris® wound gel has to be stored at room temperature between +5° and +25°C. Protect from direct sunlight.

Shelf life is 30 months. Expiry date is on the packaging. After opening the tube, the grace period is 3 months; always close tube tightly after use.

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General information – ActiMaris® Sensitiv wound rinsing solution (300ml)

1. Product profile

ActiMaris® Sensitiv wound irrigation solution is made of water with sea salt (1.2%), NaOCl (0.04%) and HOCI (0.0044%).

ActiMaris® Sensitiv wound irrigation solution in interaction with the mucous membrane, skin and wound is nascent.

Active oxygen, sodium hypochlorite and sea salt have an eliminating effect on typical sore pathogens without damaging the tissue. The decontaminating effect against MRSA / VRE tested positive. The oxidative effect of singlet oxygen (O_2) and hypertonic saline are also responsible for the neutralization of wound odors.

The basic pH-value of ActiMaris® Sensitiv wound irrigation solution is supported by the provision of OH-ions (as part of a naturally balanced redox system) in addition to the controlled release of oxygen.

Sea salt contains sodium chloride as its main ingredient (NaCl = 98%) as well as 2% other sea salts.

The increased salinity of ActiMaris® Sensitiv wound irrigation solution in turn ensures optimum hyperosmotic effect so that peri-wound swelling subsides quickly.

2. Fields of application

ActiMaris® Sensitiv wound irrigation solution is suitable:

For cleaning, humidification, decontamination, biofilm resolution, swelling and physiological debridement, to create an optimal «micro-environment» of:

- Acute mechanical wounds (skin lacerations, bites, cuts and abrasions, crack and contused) and post-operative wounds.
- Chronic wounds (pressure ulcers, arteriovenous ulcers, diabetic ulcers)
- Necrotic, malodorous wounds and ulcerating tumors, even with cavities (Hollow spaces)
- Thermal and chemical wounds (burns 1st to 3rd degree)
- entry portals of urological catheters and PEG tubes and drains
- with venous ulcers
- intrathoracic, peritoneal and ostomy rinses

For the intra-operative cleaning and rinsing of superficial wounds.

Moistening of bandages and wound dressings.

During dressing changes to loosen encrusted bandages or other encrusted wound dressings.

For the local treatment of the skin and mucosa in inflammatory and infectious processes.

For the prevention of infection at the mucosa and skin, burns and other wound types.

3. Instructions for use

ActiMaris® Forte wound rinsing solution can be applied at room temperature. In the event of strong pain, it is suggested that the necessary amount be heated to body temperature in a separate bowl before use (rinsing solution is to be heated to approximately 37°C). The duration and frequency of application should be adapted to the individual wound conditions and should so long be carried out until all deposits and necrosis are easily removed and the wound is visually clean. Depending on the wound condition, a bandage change can be done according to indications and phases, every 2 to 3 days. In case of non-infected wounds, necrosis and coatings can be mechanically removed depending on the wound condition. In cases of difficult to remove bandages it is recommended to hydrate the wound dressings with wound irrigation solution ActiMaris® (Forte or sensitivity) for at least 5 minutes until a smooth removal of bandages without trauma to the wound surface is possible. A primary surgical therapy is recommended with thick, firmly adhering fibrinous or necrotic wound deposits.

For cleaning wounds, ActiMaris® Wound Irrigation Solution soaked bandages (See also Section 8: interactions) can be combined with inert dressings or other wound dressings.

The application should preferably be in the form of the classic wet / dry phase (=moist balance). The duration of the wet phase is to be adjusted accordingly depending on the thickness and strength of the fibrin, microbial colonization, wound edge maceration or irritation.

- Stable wound environment, clean wound bed without microbial presence:
 - Recommendation: 5-10 minutes of application time
- Fibrin coating, without microbial presence, also applied to wound edge maceration or irritation:
 - Recommendation: 15 minutes of application time
- Fibrin coating, foul-smelling wounds, heavily populated by bacteria:
 - Recommendation: 20 minutes of application time

Be careful when applying the wound irrigation solution so that the bottle does not touch the skin or wound of the patient in order to prevent contamination or the spread of micro-organisms.

4. Tissue tolerance and biocompatibility

ActiMaris® Sensitiv wound rinsing solution is suitable for long-term use even on sensitive / irritated skin and mucosa due to its non-irritating and well-tolerated application. There is no inhibition of granulation and epithelialisation. Dermatologically approved.

5. Risks followed by improper use

- Improper use may cause softened wound edges
- Cross-contamination between wounds or patients with incorrect handling of the tube / bottle

6. Side effects

So far, in several years of clinical experience, the following symptoms occasionally occurred

- Slight temporary burning sensation
- Slight burning and sneezing if applied to nose area.

7. Contraindications

Since allergies can never be excluded ActiMaris® irrigation solutions or ActiMaris® gel should not be used with known or suspected allergies to any of its ingredients.

8. Restrictions on use

Pregnancy and lactation:

There is no found evidence of genotoxicity, carcinogenicity or embryotoxic effects of the ingredients. Due to the lack of studies and long-term clinical experience with pregnant and breast feeding women, ActiMaris® wound gel should be used in these cases only after an accurate medical assessment.

Infants and newborns:

Due to insufficient experience of long-term exposure, ActiMaris® wound gel should only be applied after careful medical benefit-risk assessment.

9. Interactions:

ActiMaris® Sensitiv wound rinsing solution acts as a bio-physical compound due to the content of ROS substances ($O_2 + OCl^-$) energy and oxidizing, which is why a simultaneous combination with other local, active substances and / or interactive wound dressings is in principle unnecessary. Contact with non-precious metals should be avoided. The wound can be covered with modern wound dressings after applying ActiMaris® Sensitiv wound rinsing solution.

On the other hand, a good compatibility with modern wound dressings such as silver dressings has been shown in clinical trials. Likewise, there is ample clinical evidence of good tolerability when combined with inert, non-active bandages and wound dressings such as compresses, gauze, pads, sponges, gels, hydro fibers, alginates, hydrocolloids etc. The concomitant use of ActiMaris® wound irrigation solution with the ActiMaris® Forte wound irrigation solution is possible and useful.

10. General safety instructions

ActiMaris® Sensitiv wound rinsing solution is suitable only for external wounds on skin or mucosa. Do not use for infusion or injection!

Do not swallow!

Use only untouched and undamaged containers. Keep out of reach of children!

It is possible for colored fabrics to fade after contact with ActiMaris® products.

11. Composition and technical information

Sea salt, water.

Clear, colorless, basic and hypertonic sea salt solution with characteristic smell; contains 1.2% sea salt (Sal Maris) and 0.04% NaOCl.

12. Storage specification, duration, expiry date

ActiMaris® Sensitiv wound rinsing solution has to be stored at room temperature between +5° and +25°C. Protect the bottle from direct sunlight.

Shelf life is 30 months. Expiry date is on the packaging.

The period of use is 3 months after opening the bottle. Close carefully after use.

ActiMaris AG **Manufacturer:**
 ActiMaris AG
 Technikumstrasse 14
 CH-9470 Buchs
 admin@actimaris.com
 Tel.: +41 71 505 75 25

smile medications

- Lot number
- Order number
- Follow instructions for use
- Use by: year and month
- Storage (temperature)
- Protect from direct sunlight
- Manufacturer

Date of information: 06/2017 - V5.0

Medical product: Klasse II b, CE 1250